PLACA FEMORAL DISTAL

Manual quirúrgico

Placa femoral de ángulo estable

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Distal femoral fractures often present the surgeon with a great challenge. Whereas younger patients often sustain complex fractures involving the joint with major soft tissue injury due to high-speed trauma, in geriatric traumatology the problem is one of osteoporotic fractures with reduced bone substance. Furthermore, periprosthetic femoral fractures will become increasingly important in the future. Both as regards the fracture and also because of the accompanying circumstances, implants are necessary that ensure the best possible stability and that can also be inserted minimally invasively, thus sparing the soft tissue. The distal femur plate from Königsee Implantate meets these requirements. The design was tested in an anatomical study in cadaver bone, using conventional radiology.

The plate is anatomically shaped. On the other hand, it is not applied directly to the bone so that the periosteal circulation can be preserved. Angle-stable interlocking is possible at all of the holes in the plate. By using the fluoroscopically visible guide frame, the screws can be placed in the holes in the shaft region without complex access. In the case of periprosthetic fractures with a prosthesis stem, stabilisation can be achieved by means of short screws inserted through the guide frame.

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The distal femoral plate is employed in the management of fractures of the distal femoral shaft and fractures close to or involving the joint:

- Condyle fractures with joint involvement
- Supracondylar fractures
- Distal femoral shaft fractures
- Periprosthetic femoral fractures

The plate is particularly suitable for management of 33-A fractures and 33-C fractures.

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

The postoperative treatment should be determined by the surgeon.

Primary full loading will usually not be achievable. However, stable movements should always be achieved with treatment on the motor rail. An increase to full loading is to set individually.

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In case of absent cortical support on the medial side (widespread shatter zone) in connection with expected poor bone healing the doctor has to decide about additional stabilizing intra- and post-operative measurements in order to prevent implant malfunction by long-term overload.

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
Use of planning templates assists preoperative planning. These allow the plate length to be determined before the start of the surgery. The operation is performed with the patient supine on an X-ray permeable table. If the fracture is located very far distally, it is sufficient to lower the healthy leg a little. More proximal fractures require the positioning of the healthy leg in a leg holder so that X-rays can be taken in the lateral plane also. It should be operated in order to minimise the muscle pull on the condyles, which facilitates reduction. After positioning, frontal and lateral X-rays should be taken again before the start of the operation. The entire region of the plate being implanted must be imaged.

A tourniquet can be used only in the case of smaller accesses to manage entirely distal fractures.

To treat distal femoral fractures without joint involvement, the lateral access should be selected.

The skin incision is longitudinal and the iliotibial tract is split. The vastus lateralis muscle is held forward by a Hohmann retractor. This means that the outside of the lateral femoral condyle is exposed. Intraarticular fractures with joint involvement require anterolateral access with arthrotomy. This allows reduction and retention under vision.

Apart from the very distal fractures, which allow the use of the 4-hole plate, the guide frame should be used for soft tissue sparing implantation.

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
The distal part of the guide frame fits the plate design and sits directly on the plate.

Plates and guide frames are identified with "R" and "L" for the corresponding side.

A long threaded drill bush for 3.2 mm drill is first inserted through the most distal hole of the projecting part of the guide frame and fixed in the plate.

To tighten it (caution: do not damage the guide frame), the turn pin can also be used. A retaining cap is then screwed on to ensure stability. A drill bush for Kirschner wires is then inserted through the guide attachment into the most distal hole of the part of the guide frame in contact with the plate. The plate is thus fixed to the guide frame with adequate stability.

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

In the management of fractures involving the joint, reconstruction of the articular fragment is necessary first. This can be done indirectly under fluoroscopic control or directly by opening the joint. Depending on the fracture situation, retention can be temporary with Kirschner wires or definitive with screws. The screws should be positioned so that they do not interfere with subsequent implantation of the plate.

In the metaphyseal region, reduction is facilitated by lowering the lower leg. The fracture can usually be aligned indirectly by pushing the plate proximally. In the fracture region, reduction is possible in the form of the Rep-ret-system. The threaded wire is inserted through the drill bush.

Compression can be obtained with the wing nut through the retaining cap.

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
Following dissection of the lateral femoral condyle, the plate with attached guide frame can be inserted. It is useful to prepare the way between the vastus lateralis and periosteum with the tissue dissector. The plate with guide frame can then be pushed in. It should first be guided further proximally and then back and distally until it sits correctly on the lateral condyle.

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A stab incision is now made at the level of the proximal end of the plate. The plate lengths are clearly marked on the guide frame. Depending on the length of the selected plate, a threaded tissue protection sleeve containing a trocar is inserted. This is why the guide frame holes for the most proximal plate holes of all plate lengths are threaded.


After removing the trocar, a drill bush for 1.6 mm Kirschner wires is screwed in and then stabilised with a retaining cap. A Kirschner wire is then screwed through both cortices. The proximal plate end must sit on the centre of the femoral shaft. The position should be checked by a lateral X-ray.

The system is thus closed and the plate is fixed to the femur with adequate stability.

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
For fixation of the plate it is important that it is not in a tilted position on the femur. Because of the intrinsic weight of the guide frame, the plate-guide-frame construction has a tendency to tilt away from the femoral shaft. This can be prevented by carefully pressing the distal end of the plate under vision. Through the distally inserted drill bush for Kirschner wires, a 1.6 mm wire is then drilled. This establishes the distal plate position.

Ensure by X-ray that this wire is at least 3 mm from the joint so that the screw that will be inserted later does not injure the cartilage or enter the joint. Optionally, two more Kirschner wires can also be inserted through the small holes at the distal end of the guide frame.

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The screws can now be inserted according to the preoperative planning. The screws should first be placed in the condylar region according to the part of the guide frame in contact with the plate. Selftapping angle-stable 5.5 mm cortical screws are available.

The tissue protection sleeve is inserted and the drill bush for 4.0 mm spiral drill is screwed into it.

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A preliminary hole is drilled with the 4.0 mm spiral drill. The length is measured by means of the advanced measuring aid.

Alternatively, the length can also be determined with a gauge.

To do so, the drill bush is removed and the gauge is put on directly on the plate hole.

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

If 4.5 mm non-angle-stable management is planned distally, the drill bush for 3.2 mm spiral drill must be used as well as the tissue protection sleeve. The 3.2 mm drill is used for drilling.

For angle-stable screws in general: Measurement reading = screw length.


2 mm must be added for non-angle-stable screws.

Finally, the screw is inserted with the hexagonal screwdriver. Screws are inserted the same way in all the other distal plate holes.

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For filling the shaft holes, following a stab incision the unthreaded tissue protection sleeve is advanced to the plate together with the trocar through the corresponding hole in the guide frame. After removing the trocar, the drill bush is screwed in. A hole is predrilled with the 3.2 mm drill and the attached measurement aid is used to determine screw length.

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
Alternatively, the length can also be measured with the gauge after removing the drill bush. The distance sleeve remains on the gauge. Drilling is bicortical (interior implants are exceptions). The cortical screw is screwed in through the tissue protection sleeve with the hexagonal screwdriver. Finally, the filled plate hole is identified by snapping in a blind plug.

All of the other proximal shaft holes are filled in the same way. If non-angle-stable cortical screws are used, the screws should be 2 mm longer. When angle-stable screws are used, at least 4 screws should be placed in the shaft above the fracture. If necessary, any plate holes not occupied by screws can be filled with a screw head to prevent tissue migration into the plate.

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

After all the necessary shaft screws have been inserted, the guide frame must be removed before the last screw is placed on the condyle.

The 4.0 mm drill can be used to drill through the drill guide and the screw can then be inserted.

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After placing a deep drain, the wound is closed in layers with subcutaneous drain and skin suture.

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To remove the implant after the fracture has healed, the guide frame should be used with longer plates in order to minimise soft tissue trauma.

The proximal screws can be removed through the inserted sleeves and the distal ones can be unscrewed directly.

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Instruments in extension-set, steel and titanium

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10

Disector de tejidos

Instrument for taking off the tissue

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Sleeve for femur guide frame, diameter 10.0 mm, length 149 mm, with thread

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Sleeve for femur guide frame, diameter 10.0 mm, length 149 mm, without thread

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Drill guide for K-wire D 1.6 mm, length 40 mm, ID 1.8 mm, ED 10.0 mm

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
Drill guide for K-wire D 1.6 mm, length 171 mm, ID 1.8 mm, ED 8.0 mm

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Locating sieve for instruments in extension-set

19.603.00

1

Atornillador hexagonal con mango

para tornillos D 4,5 – 7,0 mm, longitud 260 mm

Hexagonal screw-driver, with handle, for screws diameter 4.5 – 7.0 mm, length 260 mm

2.940.35/260

1

Casquillo guía de tornillos

para estribo guía femoral, D 10,0 mm

Screw holding sleeve for femur guide frame, diameter 10.0 mm

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Trocar, to transfix soft tissues, diameter 8.0 mm

10.176.12

1

Broca espiral para anclaje rápido

D 3,2 mm, longitud 305 mm

Twist drill for quick coupler, diameter 3.2 mm, length 305 mm

10.176.23

1

Broca espiral para anclaje rápido

D 4,0 mm, longitud 305 mm

Twist drill for quick coupler, diameter 4.0 mm, length 305 mm

10.176.24

1

Medidor para broca D 3,2 mm

Measure instrument, for drill diameter 3.2 mm

2.952.032

2

Medidor para broca D 4,0 mm

Measure instrument, for drill diameter 4.0 mm

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Screw gauge, with clasp, for fixed angle screws diameter 4.5 – 6.5 mm, MR 110 mm

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Wing nut for thread wire for guide frame, M 3

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D 1,6 mm, longitud 280 mm

Kirschner wire with trocar point and round end, diameter 1.6 mm, length 280 mm

6.061.16

5

OP-Set angle-stable distal femur plates

for screws 4.5/5.5 mm; large fragment; titanium; set no. 19.602.00

distal femur plate; right
Item no. 5.9482.12

hole repair stick
Item no. 9.946.05

OP-Set; Item no. 19.602.00

Item description	Item no.	Quant.
------------------	----------	--------

Implants - plates		
Angle-stable distal femur plates with combination holes in shank; right		
6 holes; length 150 mm; width 18 mm	5.9462.06	1
8 holes; length 182 mm; width 18 mm	5.9482.08	1
12 holes; length 246 mm; width 18 mm	5.9482.12	1
16 holes; length 309 mm; width 18 mm	5.9482.16	1
Angle-stable distal femur plates with combination holes in shank; left		
6 holes; length 150 mm; width 18 mm	5.9472.06	1
8 holes; length 182 mm; width 18 mm	5.9492.08	1
12 holes; length 246 mm; width 18 mm	5.9492.12	1
16 holes; length 309 mm; width 18 mm	5.9492.16	1

Instruments		
hex screwdriver with handle for screws Ø 4.5 - 7.0 mm; AF 3.5; length 215 mm	2.9406.35	1
drill guide for angle stable screwing; length 70 mm	2.977.02	2
drill guide for angle stable screwing Ø 5.5 mm; length 70 mm	2.977.15	2
drill bit for quick coupling; Ø 3.2 x 195 mm; spiral length 60 mm; double spiral	2.904.09	1
drill bit for quick coupling; Ø 4.0 x 220 mm; spiral length 60 mm; double spiral	2.904.23	1
screw forceps; self holding	2.954.01	1
gauge with clasp for screws with conical head thread Ø 4.5/6.0; measuring range 90 mm	2.953.90	1

Container		
perforated autoclavable container with inset for instruments and distal femur plates	19.601.00	1

Wire		
Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	6.031.16	5

Implants - screws		
hole repair stick		
screw head with conical thread Ø 8.0 x 5.0 mm; large fragment, titanium	9.946.05	5

Item description	Item no.	Quant.
------------------	----------	--------

Implants - screws		
Angle-stable cortical screws; Ø 4.5 mm; fully threaded; self tapping		
length 14 mm	3.545.14	6
length 18 mm	3.545.18	6
length 22 mm	3.545.22	6
length 26 mm	3.545.26	6
length 30 mm	3.545.30	6
length 34 mm	3.545.34	6
length 38 mm	3.545.38	6
length 42 mm	3.545.42	6
length 46 mm	3.545.46	6
length 50 mm	3.545.50	3
length 54 mm	3.545.54	3

Angle-stable cortical screws; Ø 5.5 mm; fully threaded; self tapping		
length 40 mm	3.549.40	2
length 45 mm	3.549.45	2
length 50 mm	3.549.50	2
length 55 mm	3.549.55	4
length 60 mm	3.549.60	4
length 65 mm	3.549.65	4
length 70 mm	3.549.70	4
length 75 mm	3.549.75	4
length 80 mm	3.549.80	4
length 85 mm	3.549.85	4
length 90 mm	3.549.90	2


Cortical screws Ø 4.5 mm; fully threaded; self tapping		
length 14 mm	3.152.14	6
length 18 mm	3.152.18	6
length 22 mm	3.152.22	6
length 26 mm	3.152.26	6
length 30 mm	3.152.30	6
length 34 mm	3.152.34	6
length 38 mm	3.152.38	6
length 42 mm	3.152.42	6
length 46 mm	3.152.46	6
length 50 mm	3.152.50	3
length 54 mm	3.152.54	3

OP-Set angle-stable distal femur plates


for screws 4.5/5.5 mm; large fragment; stainless steel; set no. 19.600.00

distal femur plate; right
Item no. 4.9482.12

hole repair stick
Item no. 8.946.05

OP-Set; Item no. 19.600.00 (shown above is the set in titanium)

Item description	Item no.	Quant.
Implants - plates		
Angle-stable distal femur plates with combination holes in shank; right		
6 holes; length 150 mm; width 18 mm	4.9462.06	1
8 holes; length 182 mm; width 18 mm	4.9482.08	1
12 holes; length 246 mm; width 18 mm	4.9482.12	1
16 holes; length 309 mm; width 18 mm	4.9482.16	1
Angle-stable distal femur plates with combination holes in shank; left		
6 holes; length 150 mm; width 18 mm	4.9472.06	1
8 holes; length 182 mm; width 18 mm	4.9492.08	1
12 holes; length 246 mm; width 18 mm	4.9492.12	1
16 holes; length 309 mm; width 18 mm	4.9492.16	1
Instruments		
hex screwdriver with handle for screws Ø 4.5 - 7.0 mm; AF 3.5; length 215 mm	2.9406.35	1
drill guide for angle stable screwing; length 70 mm	2.977.02	2
drill guide for angle stable screwing Ø 5.5 mm; length 70 mm	2.977.15	2
drill bit for quick coupling; Ø 3.2 x 195 mm; spiral length 60 mm; double spiral	2.904.09	1
drill bit for quick coupling; Ø 4.0 x 220 mm; spiral length 60 mm; double spiral	2.904.23	1
screw forceps; self holding	2.954.01	1
gauge with clasp for screws with conical head thread Ø 4.5/6.0; measuring range 90 mm	2.953.90	1
Container		
perforated autoclavable container with inset for instruments and distal femur plates	19.601.00	1
Wire		
Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	6.031.16	5
Implants - screws		
hole repair stick 	8.946.05	5
screw head with conical thread Ø 8.0 x 5.0 mm; large fragment, stainless steel	8.946.05	5

Item description	Item no.	Quant.
Implants - screws		
Angle-stable cortical screws; Ø 4.5 mm; fully threaded; self tapping 		
length 14 mm	2.545.14	6
length 18 mm	2.545.18	6
length 22 mm	2.545.22	6
length 26 mm	2.545.26	6
length 30 mm	2.545.30	6
length 34 mm	2.545.34	6
length 38 mm	2.545.38	6
length 42 mm	2.545.42	6
length 46 mm	2.545.46	6
length 50 mm	2.545.50	3
length 54 mm	2.545.54	3

Angle-stable cortical screws; Ø 5.5 mm; fully threaded; self tapping 		
length 40 mm	2.549.40	2
length 45 mm	2.549.45	2
length 50 mm	2.549.50	2
length 55 mm	2.549.55	4
length 60 mm	2.549.60	4
length 65 mm	2.549.65	4
length 70 mm	2.549.70	4
length 75 mm	2.549.75	4
length 80 mm	2.549.80	4
length 85 mm	2.549.85	4
length 90 mm	2.549.90	2

Cortical screws Ø 4.5 mm; fully threaded; self tapping 		
length 14 mm	2.152.14	6
length 18 mm	2.152.18	6
length 22 mm	2.152.22	6
length 26 mm	2.152.26	6
length 30 mm	2.152.30	6
length 34 mm	2.152.34	6
length 38 mm	2.152.38	6
length 42 mm	2.152.42	6
length 46 mm	2.152.46	6
length 50 mm	2.152.50	3
length 54 mm	2.152.54	3

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Königsee Implantate GmbH